WG-2429

[WEEE LOGO]

SERVICE.TRISTAR.EU

Tristar Europe B.V. | Jules Verneweg 87
5015 BH Tilburg | The Netherlands

EN Instruction manual

SAFETY

- By ignoring the safety instructions the manufacturer can not be held responsible for the damage.
- The appliance must be placed on a stable, level surface.
- The body fat analysis scale is not suitable for the following persons:
 - Children under 10 years and adults over 85 years
 - Persons with symptoms of fever, oedema, or osteoporosis
 - Persons undergoing treatment for dialysis
 - Persons who use a pacemaker
 - Persons who take cardiovascular medication
 - Pregnant women
 - Sportspersons who perform over 5 hours of intensive sporting activities each day
 - Persons with a heartbeat lower than 60 beats/minute

Batteries

- Do not expose the battery to high temperatures or direct sunlight. Never throw batteries into the fire. There is a danger of explosion!
- Keep batteries away from children. They are not a toy!
- Do not open the batteries by force.
- Avoid contact with metallic objects.(Rings, nails, screws etc.)there is a danger of short-circuiting!
- As a result of a short-circuit batteries may heat up considerably or even catch fire. This may result in burns.
- For your safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks, do not rub your eyes with your hands. If in contact with the liquid, the hands should be washed thoroughly and eyes rinsed out with clear water, if symptoms persist please consult a doctor.

BEFORE THE FIRST USE

- Take the appliance and accessories out of the box. Remove the stickers, protective foil or plastic from the device.
- Always weigh/use the scale on a hard and flat surface. When you are measuring, please keep still.
- It is advisable to take measurement at the same time of the day and at the same place.
- Always remove your shoes and socks, it is advisable to discard as many as possible clothing, wipe clean your feet before using and dry them well. Never use the appliance with wet feet.
- Step on the weighing platform gently and when you are measuring, please keep still.
- The maximum weight capacity is 150kg. The minimum measurable weight is 2.5kg. Graduation is 0.1kg.
- Remove battery compartment cover at the back of device and attach new CR2032 Lithium battery. Replace the cover.
- Select with the little knob aside the battery holder the kilograms, pounds or stones, on the display vershijft KG, LB or ST.

USE

- Remove the cover on the bottom of the scale and place a CR2032 lithium battery, replace the cover.
- Stand on the scale. The scale will turn on automatically.
- Always stand in the centre of the platform.
- Within seconds, personal weight is established electronically and displayed.
- Now for example you can take up your baby, if the weight is heavier than 2 kg, the scale will automatically indicate the new weight, you need to remain standing on the scale.
- Once you step off the scale after weighting, it will switch itself off automatically after a few seconds.

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CLEANING AND MAINTENANCE

- Do not immerse the appliance in water or any other liquid.
- Clean the appliance with a damp cloth. Never use harsh and abrasive cleaners, scouring pad or steel wool, which damages the device.

GARANTIE

- This product is guaranteed for 24 months granted. Your warranty is valid if the product is used in accordance to the instructions and for the purpose for which it was created. In addition, the original purchase (invoice, sales slip or receipt) is to be submitted with the date of purchase, the name of the retailer and the item number of the product. For the detailed warranty conditions, please refer to our service website: www.service.tristar.eu

ENVIRONMENT

- This appliance should not be put into the domestic garbage at the end of its durability, but must be offered at a central point for the recycling of electric and electronic domestic appliances. This symbol on the appliance, instruction manual and packaging puts your attention to this important issue. The materials used in this appliance can be recycled. By recycling of used domestic appliances you contribute an important push to the protection of our environment. Ask your local authorities for information regarding the point of recollection.

Support

You can find all available information and spare parts at service.tristar.eu!

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