EN| Instruction manual

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PT| Manual de utilizador

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SV| Bruksanvisning

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SK| Návod na použitie

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EN Instruction manual

SAFETY INSTRUCTIONS

- By ignoring the safety instructions the manufacturer can not be hold responsible for the damage.
- Please read the instruction manual carefully before using the device. Please keep these instructions, the guarantee certificate, the sales receipt and, if possible, the carton with the inner packaging.
- This appliance is not intended to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instructions concerning the use of the appliance by a person that is responsible for their safety.
- To protect children against the dangers of electrical appliances, please make sure that you never leave the appliance unattended. Therefore you have to select a storage place for the appliance where children are not able to grab it. Make sure that the cable is not hanging in a downward position.
- This appliance is only to be used for household purposes and only for the purpose it is made for.

- The appliance must be placed on a stable, level surface.
- All repairs should be made by a competent qualified repair service(*).
- Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it.
- Do not load the appliance with more than 200kg to avoid damaging the device.
- If the scales are carried or not going to be used, please store them only in a horizontal position.
- The body fat analysis scale is not suitable for the following persons:
 - Children under 10 years and adults over 85 years
 - Persons with symptoms of fever, oedema, or osteoporosis
 - Persons undergoing treatment for dialysis
 - Persons who use a pacemaker
 - Persons who take cardiovascular medication
 - Pregnant women
 - Sportspersons who perform over 5 hours of intensive sporting activities each day
 - Persons with a heartbeat lower than 60 beats/minute

Batteries

- Do not expose the battery to high temperatures or direct sunlight.Never throw batteries into the fire. There is a danger of explosion!
- Keep batteries away from children.They are not a toy!
- Do not open the batteries by force.
- Avoid contact with metallic objects.(Rings, nails, screws et cetera)there is a danger of short-circuiting!
- As a result of a short-circuit batteries may heat up considerably or even catch fire.this may result in burns.
- For your safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks, do not rub your eyes with your hands. If in contact with the liquid, the hands should be washed thoroughly and eyes rinsed out with clear water, if symptoms persist please consult a doctor.

PARTS DESCRIPTION

1. Contact surface
 2. Display
 3. On/set button
 4. Up button
 5. Down button
 6. Start button
- BEFORE THE FIRST USE**
- Take the appliance and accessories out of the box. Remove the stickers, protective foil or plastic from the device.
 - Always weigh/use the scale on a hard and flat surface. When you are measuring, please keep still.
 - It is advisable to take measurement at the same time of the day and at the same place.
 - Step on the weighing platform gently and when you are measuring, please keep still.
 - The device can only function correctly if you stand on the contact surface with bare, dry feet!
 - The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body.
 - This signal is generally safe! However, please read our special safety instructions!
 - The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data(age, gender, body weight).
 - There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values, fat is a factor that puts serious strain on the body, in particular the cardio-vascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

USE

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

Saving Personal Data

- The scale has default values that you can change, as described below.
- Carry out the following steps in sequence and within a few seconds, otherwise the device will automatically cancel the input mode.
- Press the button SET. The memory space will flash on the display (P0).
- Confirm the setting by pressing SET or select another memory space by pressing the UP/DOWN button and then SET to go to the next setting.
- The symbol for "male" or "female" will flash. Confirm the setting by pressing SET or use the UP/DOWN button to choose the symbol for your gender and then press SET to go to the next setting.
- The height setting will flash. Select your height by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The setting AGE will flash. Select your age by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing START.

Weigh on based on personal settings

- The device can only function correctly if you stand on the device with bare,dry feet. If necessary wipe the contact surface to improve the results.
- Due to the significant fluctuation of body fat the optimal testing time is between 6 p.m. And 8 p.m, since the body fat remains very consistent during this period of time.
- Press the SET button, the user number will flash (P0-P12). Use the UP/DOWN button to select your preset data.

- Wait for the user number to stop flashing.
- Stand on the contact surfaces with both feet.
- Firstly your weight will flash on the display.
- After a few seconds, the following results will be shown:

ABOUT FAT, HYDRATION AND CALORY

Fat % and hydration %

- The fat% is an index % of fat content in human body.
- The fat content in human body is of a certain level, overabundance of fat will bring bad effects on metabolism, thereby, various illnesses may invade. Please consult your doctor for more advice.
- The hydration % is an index of hydration in a human body. A to high or to low percentage is bad for your health, please advice your doctor
- Fat/hydration content fitness assessment chart
- After your weighing, you can use the following assessment charts for reference: Fig. 1

Chart age to 30 years

Fat % Female	Hydration % Female	Fat % Male	Hydration % Male	Status	Key line indication
4.0-16.0	66.0-57.8	4.0-11.0	66.0-61.2	Too lean	
16.1-20.5	57.5-54.7	11.1-15.5	61.1-58.1	Lean	
20.6-25.0	54.6-51.6	15.6-20.0	58.0-55.0	Normal	
25.1-30.5	51.5-47.8	20.1-25.5	54.9-51.9	Fat	
30.6-60.0	47.7-27.5	24.6-60.0	51.8-27.5	Total	

Chart age above 31 years

Fat %Female	Hydration % Female	Fat %Male	Hydration %Male	Status	Key line indication
4.0-20.0	66.0-55.0	4.0-15.0	66.0-58.4	Too lean	
20.1-25.0	54.9-51.6	15.1-19.5	58.3-55.3	Lean	
25.1-30.0	51.8-48.1	19.6-24.0	55.2-52.3	Normal	
30.1-35.0	48.0-44.7	24.1-28.5	52.2-49.2	Fat	
35.1-60.0	44.6-27.5	28.6-60.0	49.1-27.5	Total	

The above information is only for reference

Calorie: The data reading indicates the energy you need to sustain the daily metabolism and ordinary activities required for your weight, height, age and gender.

Error: If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.

Lo: The battery is nearly empty, replace the battery.

Err: Over load indication. The weight subject on the platform exceeds the maximum capacity of scale please step off to avoid damage.

Err-h: Low fat percentage is to low.

Err-h: High fat percentage is to high.

C: Error, please re-test to get the correct result.

CLEANING AND MAINTENANCE

- Do not immerse the appliance in water or any other liquid.
- Clean the appliance with a damp cloth. Never use harsh and abrasive cleaners, scouring pad or steel wool, which damages the device.

GUARANTEE

This product is guaranteed for 24 months granted. Your warranty is valid if the product is used in accordance to the instructions and for the purpose for which it was created. In addition, the original purchase (invoice, sales slip or receipt) is to be submitted with the date of purchase, the name of the retailer and the item number of the product.

For the detailed warranty conditions, please refer to our service website: www.service.tristar.eu

GUIDELINES FOR PROTECTION OF THE ENVIRONMENT

This appliance should not be put into the domestic garbage at the end of its durability, but must be offered at a central point for the recycling of electric and electronic domestic appliances. This symbol on the appliance, instruction manual and packaging puts your attention to this important issue. The materials used in this appliance can be recycled. By recycling of used domestic appliances you contribute an important push to the protection of our environment. Ask your local authorities for information regarding the point of recollection.

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